



## Oxygen / Hydrogen Analyzer OH-900

### General Information

The OH-900 determines oxygen and hydrogen in inorganic samples via inert gas fusion in an impulse furnace with temperatures in excess of 3,000 °C.

The OH-900 guarantees precise and fast sample analysis. The analyzer covers a wide range of applications such as metal, ceramics and other inorganic materials.

The OH-900 can be supplied with up to two infrared cells with different path lengths, accommodating both high and low level oxygen analysis. Hydrogen concentration is determined in the OH-900 by a robust and sensitive thermal conductivity cell.



### Application Examples

alloys, cast iron, ceramics, copper, refractory metals, steel, ...

### Product Advantages

- simultaneous hydrogen and oxygen determination with inert gas fusion technique
- flexible configurations and measuring ranges for O and H
- ramping and fractional analysis included
- high sensitivity IR and TC cells with low detection limits
- short analysis time
- powerful 8 kW impulse furnace for temperatures in excess of 3,000 °C
- economic analysis of grains without capsules
- easy to replace, economic upper electrode insert
- rapid, precise, accurate and reliable element determination
- powerful software (multilingual, customized display, export of results)
- single and multipoint calibration
- low maintenance
- robust design allows usage in production control and laboratory

### Features

|                       |  |
|-----------------------|--|
| Measured elements     | hydrogen, oxygen   |
| Samples               | inorganic  |
| Furnace alignment     | vertical   |
| Sample carrier        | graphite crucibles   |
| Field of application  | ceramics, engineering / electronics, steel / metallurgy                      |
| Furnace               | electrode impulse furnace (max. 8 KW), temperatures in excess of 3000 °C     |
| Detection method      | solid state infrared absorption for oxygen thermal conductivity for hydrogen |
| Typical analysis time | 120 - 180 s  |



## Oxygen / Hydrogen Analyzer OH-900

|                               |   |
|-------------------------------|---|
| Chemicals required            | magnesium perchlorate, Schuetze reagents, sodium hydroxide                              |
| Gas required                  | compressed air (4 - 6 bar / 60 - 90 psi) helium 99.995 % pure (2 - 4 bar / 30 - 60 psi) |
| Power requirements            | 3~ 400 V, 50/60 Hz, max. 8500 W   |
| <b>Dimensions (W x H x D)</b> | 55 x 80 x 60 cm   |
| Weight                        | ~ 140 kg  |
| Required equipment            | balance (resolution 0.0001g), monitor, PC   |
| Optional accessories          | carrier gas purification, gas calibration unit, voltage stabilizer 5 KVA                |

### Function Principle

#### Operation OH-900

Operation of the OH-900 is simple and safe. The samples are weighed on the interfaced balance and the weight is transferred to the linked PC. Manual weight entry is also possible. The sample is placed into the loading head and the empty graphite crucible is put on the lower electrode tip. Depending on the application the addition of auxiliary materials which lower the melting point, such as tin or nickel, may be required. The analysis time is 2 to 3 minutes, depending on the application parameters. Cell outputs are displayed in real time. All peak profiles are saved on the data base along with the results. Also all results can be transferred to a "Laboratory Information Management System" (LIMS). The OH-900 requires minimum maintenance. The particle filters and chemicals which need to be maintained are easily accessible.

#### Measuring Principle OH-900

The measuring principle of the OH-900 allows for a wide measuring range. To analyze the sample, it is weighed and placed in the sample drop mechanism. Flushing with carrier gas prevents atmospheric gas from getting into the furnace. The graphite crucible is outgassed in the impulse furnace to reduce possible contaminations (e.g. residual hydrogen). After a stabilization phase the sample is dropped into the crucible and melts. Carbon monoxide is produced by the reaction of carbon in the graphite crucible and oxygen of the sample. Nitrogen and hydrogen are released in its elemental form. The carrier gas (nitrogen) and sample gasses pass through a dust filter before entering a Schuetze reagent catalyst which converts the CO to CO<sub>2</sub>. The CO<sub>2</sub> is measured by the infrared cells to determine the oxygen content. CO<sub>2</sub> is removed chemically and the hydrogen content is measured in the thermal conductivity cell.

### Order data



## Oxygen / Hydrogen Analyzer OH-900

### ELTRA OH-900

(Please order PC, monitor, balance and consumables (starter-kit, anhydrone, sodium hydroxide, schuetze reagent) separately)

#### Measuring ranges at 1,000 mg sample weight

|            |  |
|------------|--|
| 88100-2006 | OH-900 1xO 0.1 - 200 ppm O   |
| 88100-2007 | OH-900 2xH 0.01 - 50 ppm H   20 - 1,000 ppm H  |
| 88100-2008 | OH-900 2xO 0.1 - 200 ppm O   10 ppm - 2 % O  |
| 88100-2009 | OH-900 1xO 0.1 - 200 ppm O + 2xH 0.01 - 50 ppm H   20 - 1,000 ppm H                  |
| 88100-2010 | OH-900 2xO 0.1 - 200 ppm O   10 ppm - 2 % O + 2xH 0.01 - 50 ppm H   20 - 1,000 ppm H |

#### Further measuring range combinations on request

### PC, Monitor, Balance

|            |   |
|------------|---|
| 71015      | Computer with dual core processor, 300 GB HDD, 4 GB RAM, Windows operating system, DVD-ROM, keyboard, mouse |
| 71016      | Monitor, TFT  |
| 88600-0002 | Balance (resolution 0.0001 g)   |
| 71002      | Printer   |

### Accessories

|            |   |
|------------|---|
| 22000-2001 | Gas calibration unit (integrated in analyzer)   |
| 21000      | Carrier gas purification furnace, without filling (integrated in analyzer, please order filling and quartz wool separately) |
| 72080      | Nitrogen regulator  |
| 88600-0003 | Chiller   |

### Consumables

#### Required consumables

|            |   |
|------------|---|
| 88500-0008 | Starter-kit for 500 analyses (400 graphite crucibles, 50 outer graphite crucibles, 200 inner graphite crucibles, 50 g glass wool, 50 g quartz wool) |
| 90200      | Anhydrone (magnesium perchlorate), 454 g  |
| 90210      | Sodium hydroxide, 500 g   |
| 90270      | Schuetze reagent, 100 g   |
| 90426-1001 | Filling for carrier gas purification furnace  |

#### Optional consumables

|       |                                      |
|-------|--------------------------------------|
| 90190 | Graphite crucibles, 400 pieces       |
| 90180 | Inner graphite crucibles, 100 pieces |



## Oxygen / Hydrogen Analyzer OH-900

|            |  |
|------------|--|
| 90185      | Outer graphite crucibles, 50 pieces                                  |
| 90331      | Glass wool, 454 g  |
| 90330      | Quartz wool, 50 g  |
| 91000-1001 | Calibration standard - Copper, 100 pins, 1 g each<br>~500 ppm O      |
| 91100-1001 | Calibration standard - Steel, 100 pins, 1 g each 25-40<br>ppm N      |
| 91205-1001 | Calibration standard - Titanium, 100 pins, 0.1 g each<br>10-35 ppm H |
| 91400-1001 | Calibration standard - Steel, 100 pins, 1 g each 0.5 - 1<br>ppm H    |
| 92610      | Tube of high vacuum grease   |
| 90870      | Cooling agent, 0.5 l   |

### Spare and Wear Parts

|       |                                   |
|-------|-----------------------------------|
| 31250 | Upper electrode                   |
| 31360 | Graphite tip                      |
| 71010 | Brush                             |
| 71035 | Cleaning brush / furnace brush    |
| 11062 | Reagent tubes 160x16 mm, 2 pieces |
| 11064 | Reagent tubes 280x16 mm, 2 pieces |
| 20040 | Catalyst tube                     |
| 31365 | Graphite tip holder               |
| 77033 | Circuit breaker 32 A              |